



INFRARED THERMOMETER AIR

**Infrared Non-contact Body Thermometer
for Personal Protection, Home, School,
Workplace, Public Places or Events**

Features

- High-Precision infrared sensors
- Precise non-contact rapid measurements
- Easily switch between Celsius (°C) or Fahrenheit (°F)
- Switch between body or surface/object temperature modes
- Automatically saves the last 32 temperature readings
- Clear LCD digital display with back light
- Auto power-off

Usage Instructions

- Insert two AAA batteries (not included).
- Press SET and hold for 3 seconds to switch between object and body temperature mode.
- Point at surface. For best measurements in body temperature mode, point at the middle of the forehead or behind the ear lobe and keep the thermometer vertically about 1-2" away.
- Press the measurement button (pull trigger) until temperature shows on the LCD display. Readings take approximately 1 second and are confirmed with an audible "tick" sound.

Technical Specifications

Measuring method: Forehead
Measuring distance: 1.18"-2" (3-5cm)
Measuring tolerance:
89.6°F-94.8°F ($\pm 32.5^\circ\text{F}$), 95°F-107.6°F ($\pm 32.3^\circ\text{F}$), 107.8°F-108.3°F ($\pm 32.5^\circ\text{F}$)
32°C-34.9°C ($\pm 0.3^\circ\text{C}$), 35°C-42°C ($\pm 0.2^\circ\text{C}$), 42.1°C-42.4°C ($\pm 0.3^\circ\text{C}$)
Display resolution: 32.18°F (0.1°C)
Operating environment: 50°F~104°F (10°C~40°C) Humidity \leq 85%
Transportation storage environment: 32°F~122°F (0°C~50°C) Humidity \leq 85%
Power: DC3V (2x AAA Batteries)
Battery capacity reminder : Low battery alert
Backlight: Highlight backlight
Unit of measurement: Celsius (°C) / Fahrenheit (°F)
Automatic shut-down: 15 seconds
Size: 5.8" X 3" X 1.57" (149mm X 77mm X 40mm)
Weight: 0.24 lb (109g)

Maintenance Instructions

- Please keep the inside of the sensor and probe clean, otherwise it will affect the measurement accuracy.
- Surface cleaning: Wipe surface dirt with a clean soft cloth or cotton swab with a little rubbing alcohol or water.
 - Sensor cleaning: Use a clean soft cloth or cotton swab with a little rubbing alcohol to gently wipe the inside of the probe or the top of the sensor.
 - Please read the instructions carefully before use.
 - Do not immerse the thermometer in any liquid, and do not leave in extreme temperatures for a prolonged period.
 - Dropping and/or mixing with sharp objects is not recommended.
 - Do not use in environments with strong electromagnetic interference.
 - Keep the thermometer out of the reach of children.
 - It is recommended to practice several times to become familiar with the measuring method, and try not to change the default factory settings.
 - Measurement results are not a substitute for physician diagnosis.